

LAMINECTOMY PUNCH KÜLLING.

**UP 40°, DETACHABLE.**



## LAMINECTOMY PUNCH KÜLLING.

**UP 40°, DETACHABLE.**

### **Use and benefit**

Laminectomy punches are primarily used in spinal surgery for the removal of small portions of soft tissue, cartilage and bone.


**Laminectomy Punch Külling** – The innovative concept and ergonomic design of the handle provides a much more precise control of the tip, as well as a smooth, effort-saving and uniform movement during the punching process.

Theoretical and empirical studies prove, what can be felt by simply comparing the Külling with a conventional Kerrison laminectomy punch. The new handle enables a smooth, uniform movement, thus enhancing tip control, as well as reducing the effort required.


**Learn more** [www.ulrich-swiss.ch](http://www.ulrich-swiss.ch)

The **laminectomy punch Külling** is available in the following models:


**Laminectomy punch Külling, up 40°, detachable, shank length 200 mm**



Laminectomy punch Külling, up 40°, detachable, shank length 230 mm



Laminectomy punch Külling, up 40°, detachable, shank length 270 mm



Laminectomy punch Külling, up 40°, detachable, shank length 300 mm



length  
300 mm



flat foot U07-131-..



2mm 3mm 4mm

- |      |            |
|------|------------|
| 2 mm | U07-131-20 |
| 3 mm | U07-131-30 |
| 4 mm | U07-131-40 |



standard foot U07-130-..



1mm 2mm 3mm 4mm 5mm

- |      |            |
|------|------------|
| 1 mm | U07-130-10 |
| 2 mm | U07-130-20 |
| 3 mm | U07-130-30 |
| 4 mm | U07-130-40 |
| 5 mm | U07-130-50 |

The BlackLine.

Laminectomy punch Külling, BlackLine, up 40°, detachable, shank length 200 mm




length  
200 mm






flat foot U07-1219-..

	2 mm	U07-1219-20
	3 mm	U07-1219-30
	4 mm	U07-1219-40

2mm 3mm 4mm




standard foot U07-1209-..

	1 mm	U07-1209-10
	2 mm	U07-1209-20
	3 mm	U07-1209-30
	4 mm	U07-1209-40
	5 mm	U07-1209-50

1mm 2mm 3mm 4mm 5mm

The BlackLine.

Laminectomy punch Külling, BlackLine, up 40°, detachable, shank length 230 mm





The BlackLine.

Laminectomy punch Külling, BlackLine, up 40°, detachable, shank length 270 mm



The BlackLine.

Laminectomy punch Külling, BlackLine, up 40°, detachable, shank length 300 mm



Additional models on request.

**ULRICH AG.**  
**ST. GALLEN.**

MÖVENSTRASSE 12  
CH-9015 ST. GALLEN

T +41 71 314 62 62  
F +41 71 314 62 99

[www.ulrich-swiss.ch](http://www.ulrich-swiss.ch)  
[info@ulrich-swiss.ch](mailto:info@ulrich-swiss.ch)

